

INDIAN INSTITUTE OF TECHNOLOGY KANPUR
DEPARTMENT OF CHEMICAL ENGINEERING
PG Research lab
STANDARD OPERATING PROCEDURE
FE-SEM

1. Startup Procedure (from STANDBY condition)

- 1.1 Open Nitrogen Cylinder valve, ensure the regulator output pressure is $> 5 \text{ kg/cm}^2$
- 1.2 Press STANDBY button *to reinitiate*
- 1.3 Press PUMP *to vacuum the chamber*
- 1.4 Wait till
Column Pressure comes below 2.0 E-03 Pa
Chamber pressure comes below 2.0 E-03 Pa
The pressure indication bars become 'green'
- 1.5 Press VENT to vent the chamber. (Chamber for sample insertion opens only after completion of venting process i.e. chamber pressure comes to atmospheric)

NOTE: Make the sample ready before venting the chamber so that the chamber remains open for least period of time.
- 1.6 After fixing sample properly (hand tight the screws) close the chamber and press PUMP button holding the chamber door inward.
- 1.7 Release the chamber door as pressure reaches $< 1.0 \text{ E +04}$
- 1.8 Wait for pressure reaches desirable level as stated in 1.4
- 1.9 Wait 5 mins
- 1.10 Press BEAM ON
- 1.11 *Now the instrument is ready for sample measurement*

2. **Standby procedure**

- 2.1 Save all images taken for measurement
- 2.2 Set HV to 1.00 KV
- 2.3 Set Beam intensity to 10.00
- 2.4 Set Speed to 2-3 level
- 2.5 Set Magnification to 1.00 kx
- 2.6 After stabilization to HV= 1.00 KV press BEAM ON button *to off the beam*
- 2.7 Press STANDBY *to standby the system*
- 2.8 Close Nitrogen Cylinder Valve

Emission Parameter Adjustment

1. Log off
2. Log in as 'service'
3. Password – ivaiva
4. From menu > SEM > FEG HV Control
5. Set Limits

6. Set Maximum Emission
7. Press Enter

During restart from shut off

1. After Gun Pressure – Green
2. Service login
3. From Menu > SEM > FEG HV Control Panel
4. Power On

Other Troubleshooting

If beam becomes off immediately and GUN Pressure is RED check cylinder pressure

Any error shown in health status,
> Help > Send Log file > File – save and send to
dkpanda@advscientific.com